



## Hg5<sup>®</sup> Series of Amalgam Separators

We've got amalgam separation down to a science. Removing mercury from the dental practice (and the environment) is a good thing - and we're way ahead of the game. When it comes to Amalgam Separation, no one knows more than Solmetex. 21 years ago, our experience in water science led us to create a system that easily and affordably removes amalgam waste from dental practices. Our signature product the Hg5 is the industry standard and has won more awards than all competitors combined. As a 'total solution provider,' Solmetex takes care of waste handling and recycling so you don't have to. Good for you and good for the environment. Solmetex – clearly the leader in amalgam separation

Solmetex provides the complete solution:

- ISO 11143 certified at greater than 99%
- Professionally engineered, robust, simple design.
- Easy to operate and maintenance free.
- Flexible design for hard-to-fit applications.
- Functional for wet or dry vacuum systems (install before pump on wet vacuum system, install before tank on dry vacuum system).

## Product Awards



dental townie choice awards.2015



Amalgam Containment System

1.800.216.5505  
www.solmetex.com

## Hg5HV-001

### HG5 High Volume Amalgam Separator

#### Specification High Volume System:

- 18.5”h x 21”w x 22.5”d.
- 3” inlet and outlet piping.
- Wall mount bracket.
- Services 11 - 20 chairs.

## Hg5-001

### Hg5 Amalgam Separator

#### Specification Standard:

- 29”h x 10”w x 8”d.
- 1-1/2” inlet and outlet piping.
- Fernco® couplers.
- Floor / Wall mount bracket.
- Services 1 - 10 chairs.

## Hg5-Mini

### Hg5 Mini Amalam Separator

#### Specification Mini System:

- 21.5”h x 11”w x 8”d.
- 3/4” inlet and outlet piping.
- Fernco® couplers.
- Floor / Wall mount bracket.
- Services 1 - 4 chairs.

## Hg5-002CR

### Collection Container with Recycle Kit

#### Complete container & recycle kit:

- Collection container
- Proper packaging for shipping
- Recycling
- Recycling documentation
- One size fits all systems



To order your Solmetex Hg5 System, please contact your local dental product dealer. If you don't have a dealer, contact SolmeteX for information at 1-800-216-5505 or visit <http://www.solmetex.com> Click on "dealer locator" to find the dental product dealer closest to you.

#### Easy

- Maintenance free
- Collection container changed with no tools
- Two-minute container change out
- Clear collection container lets you see when full

#### Flexible

- Free standing or wall mount
- Wet or dry vacuum systems
- Multiple sizes to fit your practice

#### Recycle

- Ship container via UPS to our contracted EPA certified mercury facility
- Recycle Kit includes container and recycle packaging plus pre-paid shipping Certificate provided on-line

#### Support

- Expert, responsive, service oriented



Solmetex® recommends Sultan PUREVAC® SC line Cleaner  
Formulated for systems with amalgam separators  
Manufactured by Sultan Healthcare - Available through your dental dealer

