


## We have amalgam separation down to a science.

The Best Just Got Better...Solmetex is proud to introduce the NXT Hg5 Amalgam Separator. The Hg5 series of Amalgam Separators add value to the equipment room by extending the life of your vacuum system and is virtually maintenance free.

### Solutions

The NXT Hg5 provides the complete solution:

- New compact design, internal manifold, ideal for hard to fit areas
- Functional with all wet and dry vacuums
- Environmentally friendly packaging
  - ISO 11143 certified by NSF International



#### Easy

- Clear system, easy to inspect
- Two-minute container change out
- No tools necessary

#### Support

- Compliance center
- Online certificates
- Auto-ship program

#### Flexible

- Versatile, wall or floor mounted
- Compatible with wet or dry vacuums
- 3 different sizes to fit your practice

# Ordering Information

To order your Solmetex NXT Hg5 System, please contact your local dental product dealer. If you don't have a dealer, visit [www.solmetex.com](http://www.solmetex.com) Click on "dealer locator" to find the dental product dealer closest to you or contact Solmetex for information at **1-800-216-5505**.



**NXT Hg5 Mini  
Amalgam Separator**

- Item # **NXT-Hg5-Mini**  
Specifications:
- 25" h x 11" w x 8" d
  - 3/4" inlet and outlet piping
  - Fernco® couplers
  - Floor / Wall mount bracket
  - Services 1 - 4 chairs



**NXT Hg5  
Amalgam Separator**

- Item # **NXT-Hg5-001**  
Specifications:
- 30" h x 11" w x 8" d
  - 1-1/2" inlet and outlet piping
  - Fernco® couplers
  - Floor / Wall mount bracket
  - Services 1 - 10 chairs



**Hg5 High Volume  
Amalgam Separator**

- Item # **Hg5HV-001**  
Specifications:
- 28" h x 18" w x 22.5" d
  - 3" inlet and outlet piping
  - Wall mount bracket
  - Services 11 - 20 chairs



**NXT Hg5  
Collection Container**

Item #  
**NXT-Hg5CAN-002**



**Hg5  
Collection Container**

[legacy] Item #  
**Hg5-002**



**RAMVAC Utility  
Collection Container**

Item #  
**Hg5-002-RVC**

**Solmetex – clearly the leader in amalgam separation.**

The Hg5 amalgam separator is the industry standard and has won more awards than all competitors combined.

