



Vacuum Line Cleaning Solution pH and Hg5® Performance November 2017

The EPA “Dental Rule” on the proper handling of amalgam waste along with the adoption of the ADA’s Best Management Practices (BMPs) went into effect nationwide on July 14, 2017. The rule states that only non-chlorine or non-oxidizing disinfectants and neutral vacuum line cleaners with a **pH >6 and <8** are to be used.

Please check back later as we update our list of recommended vacuum line cleaners fitting the new EPA regulation.

~Team Solmetex

