

We've got amalgam separation down to a science.

- Experience in water technology
- Created a system that is affordable
- Removing mercury from the practice
- Signature product – Hg5 Series
- Setting the industry standard
- Total solution provider



Solutions.

Solmetex provides the complete solution:

- ISO 11143:2008 certified by NSF International
- Professionally engineered, robust, simple design.
- Easy to operate and easy to use.
- Flexible design for hard-to-fit area.
- Functional for wet or dry vacuum systems.

Solmetex – clearly the leader in amalgam separation

The Hg5 amalgam separator is the industry standard and has won more awards than all competitors combined.



Amalgam Containment System





Hg5 Mini



Hg5



Hg5 HV



Hg5 Collection Container with recycling Kit

Ordering information

Hg5 Mini Amalgam Separator

Item #: Hg5-Mini

Specification Mini System:

- 21.5”h x 11”w x 8”d.
- 3/4” inlet and outlet piping.
- Fernco® couplers.
- Floor / Wall mount bracket.
- Services 1 - 4 chairs.

Hg5 Amalgam Separator

Item #: Hg5-001

Specification Standard:

- 29”h x 11”w x 8”d.
- 1-1/2” inlet and outlet piping.
- Fernco® couplers.
- Floor / Wall mount bracket.
- Services 1 - 10 chairs.

HG5 High Volume Amalgam Separator

Item #: Hg5HV-001

Specification High Volume System:

- 28”h x 18”w x 22.5”d.
- 3” inlet and outlet piping.
- Wall mount bracket.
- Services 11 - 20 chairs.

Collection Container with Recycle Kit

Item #: Hg5-002CR

Complete container & recycle kit:

- Collection container
- Proper packaging for shipping
- Recycling
- Recycling documentation
- One size fits all systems

Easy

Clear system, easy to inspect
Two-minute container change out
No tools necessary

Flexible

Versatile, wall or floor mounted
Compatible with wet or dry vacuums
3 different sizes to fit your practice

Recycle

Ship container via UPS
EPA certified mercury facility
Recycle Kit includes container and recycle packaging

Support

Compliance center
Online certificates
Auto-ship program

To order your Solmetex Hg5 System, please contact your local dental product dealer. If you don't have a dealer, visit <http://www.solmetex.com> Click on “dealer locator” to find the dental product dealer closest to you or contact Solmetex for information at 1-800-216-5505.